

# Fact Sheet

## \_dataMaintenance.BRAIN

**BIZERBA**

With the master data maintenance tool\_dataMaintenance.BRAIN you centrally create and maintain article data, incl. parameters for labels, and send them to Bizerba devices. By means of the quality gate you ensure that only controlled data are sent to production.



**More Information**  
High performance  
with a wide range of  
solutions and options.

### System requirements

<b>Computer and processor</b>	Processor with 1000 MHz or above recommended i5, i7 or Xeon	<b>Supported databases</b>	Microsoft SQL Server: <ul style="list-style-type: none"> <li>• SQL Server 2008</li> <li>• SQL Server 2008 Express*</li> <li>• SQL Server 2008 R2</li> <li>• SQL Server 2008 R2 Express*</li> <li>• SQL Server 2012</li> </ul>
<b>Display</b>	Monitor with a resolution of 1024 x 768 or above		
<b>Main storage</b>	About 300 MB		
<b>Supported operating systems</b>	Microsoft Windows® (x32/x64) <ul style="list-style-type: none"> <li>• Server 2008</li> <li>• Server 2008 R2</li> <li>• Server 2012</li> <li>• Server 2012 / R2</li> <li>• Windows Vista</li> <li>• Windows 7 <ul style="list-style-type: none"> <li>- Professional</li> <li>- Ultimate</li> <li>- Enterprise</li> </ul> </li> <li>• Windows 8* / 8.1*</li> <li>• Windows 10 <ul style="list-style-type: none"> <li>- Professional</li> <li>- Enterprise</li> </ul> </li> </ul> (* Windows 8 RT is not supported.)		(* Limitations of Express version due to database limitations (1 CPU, 1 GB RAM, database 2, or max. 4 GB)
<b>Drive</b>	DVD drive	<b>Article administration</b>	Text administration, also statistic texts, barcodes, device data, automate parameter
<b>Edition Basic</b>	up to 9 individual connections	<b>Import export</b>	CSV import und export of article data from _dataMaintenance.BRAIN 2.20
<b>Edition Standard</b>	up to 39 individual connections	<b>Supported Bizerba devices</b>	<ul style="list-style-type: none"> <li>• Labeler GX devices from device version 8.60</li> <li>• Checkweigher CX devices from device version 8.60</li> <li>• Inspection systems BVS-I, BVS-L</li> </ul>
<b>Edition Professional</b>	no limitation of individual connections	<b>Definition of barcode contents</b>	QR Code, Datamatrix code, among others
		<b>Nutrition value charts</b>	Input of nutrition value data per 100g for every article, GDA and nutrition values of the portion size are calculated automatically.

### Members of the Bizerba Group

Bizerba GmbH & Co. KG  
Wilhelm-Kraut-Straße 65  
72336 Balingen  
Deutschland  
T +49 7433 12-0  
F +49 7433 12-2696  
marketing@bizerba.com  
www.bizerba.com